



Ansell**CARES**<sup>SM</sup>  
Education. Evidence. Engagement.

## **RADIATION EXPOSURE**

---

### **The Consequences Faced By Healthcare Professionals**

**CLINICAL BULLETIN**




# A CLOSE LOOK AT HOW THE C-ARM FLUOROSCOPY GENERATES RADIATION

Healthcare workers are exposed to greater amounts of ionizing radiation due to caring for patients during their treatments. Exposure may come from Radiology, mammography, computed tomography (CT), and fluoroscopy.




Radiology




Computerized Tomography (CT)



Fluoroscopy



Mammography



- A X-ray tube:** X-ray tube, the beam source, emits X-rays that penetrate the patient.
- B Primary radiation beam:** emitted from the beam source, X-ray tube.
- C Scattered radiation field:** secondary ionizing radiation spreads out in different directions when the radiation beam is intercepted by objects, such as a patient.
- D Image detector:** Captures X-ray and converts to an image displayed on a monitor.

- ✓ Primary radiation is a major source of dose to patients
- ✓ Scattered radiation is a major source of exposure to healthcare professionals
  - Reduction of scattered radiation is the main aim of RAGs

## Biological Effects of Exposure to Ionizing Radiation to Both Patients and Healthcare Professionals

### Deterministic Effects<sup>5</sup>

- ✓ Occur only **after a threshold dose** has been surpassed before a biological response is seen
- ✓ Severity is dose dependent above the threshold

#### VISIBLE



Cataracts



Skin reddening (erythema)



Skin decay (necrosis)

#### NON-VISIBLE



Temporary or permanent sterility

### Stochastic Effects<sup>6</sup>

- ✓ Occur by statistical chance and **not dependent on threshold dose:**
- ✓ No radiation dose is considered to be 'safe' for radiation protection purposes\*

#### VISIBLE



Skin cancer



Lung cancer

#### NON-VISIBLE



Genetic effects such as mutations that can be transmitted to offspring

\*In accordance with the principle of ALARA



## Risks of Exposure to Radiation are Cumulative, Long Term and Not Immediately Visible



A retrospective study revealed a **29% incidence of cancer among ortho surgeons exposed to medical radiation** compared to 4% among their unexposed peers<sup>7</sup>

- ✓ Typical lag period between radiation exposure and cancer diagnosis is at least 5 years. In most cases, **may be 1 or 2 decades** or longer<sup>2</sup>
- ✓ Latent period for radiogenic skin cancer ranges from **2 to 65 years** after exposure, with an estimated median latency of 20 to 45 years<sup>3</sup>

## General Principles of Radiation Protection and What It Means to Healthcare Professionals

### ALARA is the Most Referenced Standard of Practice

#### ALARA

The principle of **As Low As Reasonably Achievable (ALARA)**<sup>8</sup> is used by key safety and occupational agencies to minimize radiation does by employing all reasonable methods.

- Central to many radiation safety programs

- 1 Less TIME:** The shorter the time spent near a radiation source, the lower the amount of radiation exposure received.
- 2 More DISTANCE:** The closer to the radiation source, the greater the chances of damaging the body: When the working distance from a radiation source is increased by a factor of 2, the dose received from the source is reduced by a factor of 4.
- 3 More SHIELDING:** Increasing the shielding around a radiation source decreases the exposure.

## WHERE ARE HANDS IN SCATTERED RADIATION FIELD?




### SPOT THE PROXIMITY OF THE HANDS!

Hands are often close to the scattered radiation field in image guided interventional procedures due to the need to position patients.


#### O.R. PROCEDURES

Examples: Spinal surgery and Orthopedic surgery




#### CATH LAB/ INTERVENTIONAL RADIOLOGY

Examples: Angiograms, Catheter insertion, Placement of devices and Barium X-rays





**Hands are often the most exposed, yet least protected to scattered radiation?**



**IONIZING RADIATION ATTENUATION**

**Radiation Attenuation**

- Attenuation is always relative to the radiation detected with no material between the source and the detector
- It is a measure how much the radiation interacts (is lost) as it goes through a material
- Radiation can be “lost” between source and detector:
  - » When it is absorbed (energy deposited in the material)
  - » When it is scattered in a different direction, so it does not reach the detector



## KEY TAKEAWAYS



**Hands are the most exposed, yet least protected:** often closest to the primary radiation beam and scattered radiation field.

- **Hands have 13 times greater exposure** than the sum total of all other protected body parts<sup>9</sup>
- Body parts (chest, thyroid), with much less exposure, are protected with protective shields but not hands



Even though healthcare professionals may not exceed the annual permissible rate (as guided by International Commission of Radiation Protection) for hand dose, **the long term, cumulative impact of hand exposure should not be dismissed.**

- **Even the smallest doses could put a user at risk** of cancer and genetic mutations
- Dose limits **DO NOT** and **CANNOT** define a demarcation between 'safe' and 'unsafe'



**Radiation attenuation gloves are not widely used:** consistent and constant education and awareness of the level of hand dose exposure compared to other body parts and the ensuing stochastic risks are needed to **shift behavior from 'rarely or never use' to 'use regularly'.**



There are no laws mandating the use of protective shields but **ALARA is the principle of best practice** most referred to by protective agencies and healthcare occupational bodies.

### References:

1. Choudhary S. Deterministic and Stochastic Effects of Radiation. *Canc Therapy & Oncol Int J*. 2018; 12(2): 555834.
2. Lin EC. Radiation risk from medical imaging. *Mayo Clin Proc*. 2010;85(12):1142-1146.
3. Meibodi NT, Maleki M, Javidi Z, Nahidi Y. Clinicopathological evaluation of radiation induced basal cell carcinoma. *Indian J Dermatol*. 2008;53(3):137-139.
4. Fritzsche H, Phillips M. Electromagnetic radiation. *Encyclopedia Britannica*. Updated July 23,2020. Accessed September 18, 2023. Website. <https://www.britannica.com/science/electromagnetic-radiation>.
5. Coeytaux K, Bey E, Christensen D, Glassman ES, Murdock B, Doucet C. Reported radiation overexposure accidents worldwide, 1980-2013: a systematic review. *PLoS One*. 2015;10(3):e0118709. Published 2015 Mar 19.
6. Kaplan DJ, Patel JN, Liporace FA, Yoon RS. Intraoperative radiation safety in orthopaedics: a review of the ALARA (As low as reasonably achievable) principle. *Patient Saf Surg*. 2016;10:27. Published 2016 Dec 12.
7. Bratschitsch G, Leitner L, Stücklschweiger G, et al. Radiation Exposure of Patient and Operating Room Personnel by Fluoroscopy and Navigation during Spinal Surgery. *Sci Rep*. 2019;9(1):17652. Published 2019 Nov 27.
8. ICRP, 2007. The 2007 Recommendations of the International Commission on Radiological Protection. ICRP Publication 103. Ann. ICRP 37 (2-4). Updated 2007. September 18 2023. ICRP Website. <https://www.icrp.org/publication.asp?id=ICRP%20Publication%20103>.
9. Hoffer CE, Ilyas AM. Fluoroscopic radiation exposure: are we protecting ourselves adequately?. *J Bone Joint Surg Am*. 2015;97(9):721-725.
10. Wagner LK, Mulhern OR. Radiation-attenuating surgical gloves: effects of scatter and secondary electron production. *Radiology*. 1996;200(1):45-48.

➤ **For more information or additional clinical resources, please visit: [www.ansell.com/AnsellCARES](http://www.ansell.com/AnsellCARES)**

Ansell, ® and ™ are owned by Ansell Limited or one of its affiliates. © 2023 Ansell Limited. All rights reserved.