Bahaya Sarung Tangan Berbedak

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EMEA

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium

Australia

Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
Australia

Malaysia

Ansell Global Trading Center Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6,
63000 Cyberjaya, Malaysia

Referensi

-
- i Edlich RF, Long WB 3rd, Gubler DK, et al. Dangers of cornstarch powder on medical gloves: seeking a solution. *Ann Plast Surg.* 2009;63(1):111-115.
- ii Baker CM. A system in need of repair – medical device regulation: the example of latex medical gloves. *Synesis: A Journal of Science, Technology, Ethics, and Policy.* 2013;4:G32-G39. http://www.synesisjournal.com/vol4_g/Baker_2013_G32-39.pdf. Diakses 6 November 2015.
- iii Gibilisco PA. Clinical perils: latex allergy. <http://www.icumed.com/media/84928/M1-1252%20Clinical%20Perils%20of%20Latex%20Rev.03PQ.pdf>. Diakses 6 November 2015.
- iv Tarlo SM, Sussman G, Contala A, Swanson MC. Control of airborne latex by use of powder-free latex gloves. *J Allergy Clin Immunol.* 1994;93(6):985-989.
- v American Latex Allergy Association. Statistics. <http://latexallergyresources.org/statistics>. Diakses 6 November 2015.
- vi Gibilisco PA. Clinical perils: latex allergy. <http://www.icumed.com/media/84928/M1-1252%20Clinical%20Perils%20of%20Latex%20Rev.03PQ.pdf>. Diakses 6 November 2015.
- vii Gibilisco PA. Clinical perils: latex allergy. <http://www.icumed.com/media/84928/M1-1252%20Clinical%20Perils%20of%20Latex%20Rev.03PQ.pdf>. Diakses 6 November 2015.
- viii US Food and Drug Administration. Medical glove powder report. September 1997. <http://www.fda.gov/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm113316.htm>. Diakses 6 November 2015.
- ix Newsom SW, Shaw P. Airborne particles from latex gloves in the hospital environment. *Eur J of Surgery Suppl.* 1997;(579):31-33.

United States

Ansell Protective Products Inc.
111 Wood Avenue South
Suite 210
Iselin, NJ 08830, USA

EMEA

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium

Australia

Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
Australia

Malaysia

Ansell Global Trading Center Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6,
63000 Cyberjaya, Malaysia

- x Public Citizen. Petition to FDA to ban powdered and latex surgeon's and patient examination gloves. 25 April, 2011. <http://www.citizen.org/documents/1945.pdf>. Diakses 6 November 2015.
- xi European Commission. Guideline on medical devices. Implications of the medical devices directives in relation to medical devices containing natural rubber latex: a guide for manufacturers and notified bodies. http://ec.europa.eu/health/medical-devices/files/meddev/2_5_9rev_latex_en.pdf. Diakses 23 November 2015.
- xii US Food and Drug Administration. Medical glove powder report. <http://www.fda.gov/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm113316.htm>. Diakses 23 November 2015.
- xiii US Food and Drug Administration. Medical glove powder report. <http://www.fda.gov/medicaldevices/deviceregulationandguidance/guidancedocuments/ucm113316.htm>. Diakses 23 November 2015.
- xiv Wolfe SM, Sullivan T. Petition to ban cornstarch powder on latex gloves. http://www.citizen.org/publications/print_release.cfm?ID=6629. Diakses 23 November 2015.
- xv Edelstam G, Arvanian L, Karlsson G. Glove powder in the hospital environment – consequences for healthcare workers. *Int Arch Occup Environ Health*. 2002;75(4): 267–271.
- xvi Giercksky K. Misdiagnosis of cancer due to multiple glove powder granulomas. *Eur J Surg Suppl*. 1997; 163(579): 11-14.
- xvii
https://www.federalregister.gov/documents/2016/12/19/2016-30382/banned-devices-powdered-surgeons-gloves-powdered-patient-examination-gloves-and-absorbable-powder?source=govdelivery&utm_medium=email&utm_source=govdelivery.
Diakses 9 Januari 2016.
- xviii
https://www.federalregister.gov/documents/2016/12/19/2016-30382/banned-devices-powdered-surgeons-gloves-powdered-patient-examination-gloves-and-absorbable-powder?source=govdelivery&utm_medium=email&utm_source=govdelivery.
Diakses 9 Januari 2016.

United States

Ansell Protective Products Inc.
111 Wood Avenue South
Suite 210
Iselin, NJ 08830, USA

EMEA

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium

Australia

Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
Australia

Malaysia

Ansell Global Trading Center Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6,
63000 Cyberjaya, Malaysia

- xix Allmers H, Schmengler J, John SM, Skudlik C. The disappearance of NRL-allergy in Germany and Europe. <http://www.mrepc.com/cn/publications/papers/NRL-Allergy%20in%20Germany%20and%20Europe.pdf>. Diakses 6 November 2015.
- xx Vandenplas O, Larbanois A, Vanassche F, et al. Latex-induced occupational asthma: time trend in incidence and relationship with hospital glove policies. *Allergy*. 2009;64(3):415-420.
- xxi LaMontagne AD, Radi S, Elder DS, Abramson MJ, Sim M. Primary prevention of latex related sensitisation and occupational asthma: a systematic review. *Occup Environ Med*. 2006;63(5):359-364.
- xxii Edlich RF, Long WB, Gubler DK, et al. Dangers of cornstarch powder on medical gloves - seeking a solution. *Ann Plast Surg*. 2009; 63(1):822-826.
- xxiii Truscott W. The citizen's petition to ban cornstarch powder on medical gloves. In: Edlich RF, Garrison JA, Smith HN, eds. *Deadly Powder on Medical Gloves: A Wake-Up Call to the Food and Drug Administration*. Bloomington, IN: iUniverse; 2012: 48.
<https://books.google.com/books?id=9zM71zEaH3QC&pg=PA48&lpg=PA48&dq=countries+that+have+banned+the+use+of+powdered+medical+gloves&source=bl&ots=7OpI4XgUZU&sig=XJy6uPtFHOWb9jLUtGXUatWlxtY&hl=en&sa=X&ved=0CDsQ6AEwBGoVChMI87XlkbWDyQIVwXg-Ch3cjQ3u#v=onepage&q=countries%20that%20have%20banned%20the%20use%20of%20powdered%20medical%20gloves&f=false>. Diakses 9 November 2015.
- xxiv Safe Work Australia. Factors contributing to the development of occupational contact dermatitis and occupational contact urticaria. Mei 2014.
<http://webcache.googleusercontent.com/search?q=cache:7UQX9wT4x6YJ:www.safeworkaustralia.gov.au/sites/SWA/about/Publications/Documents/851/Factors-contributing-contact-dermatitis.docx+&cd=5&hl=en&ct=clnk&gl=us>. Diakses 9 November 2015.
- xxv Edlich RF, Long WB, Gubler DK, et al. Dangers of cornstarch powder on medical gloves - seeking a solution. *Ann Plast Surg*. 2009; 63(1):822-826.

United States

Ansell Protective Products Inc.
111 Wood Avenue South
Suite 210
Iselin, NJ 08830, USA

EMEA

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium

Australia

Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
Australia

Malaysia

Ansell Global Trading Center Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6,
63000 Cyberjaya, Malaysia

xxvi Mayo has had little problems with latex allergies.

http://www.postbulletin.com/mayo-has-had-little-problems-with-latex-allergies/article_5d9046e2-55ce-57dc-b427-fe01a22680_e1.html. Diakses 6 November 2015.

xxvii Hunt LW, Kelkar P, Reed CE, Yunginger JW. Management of occupational allergy to natural rubber latex in a medical center: the importance of quantitative latex allergen measurement and objective follow-up. *J Allergy Clin Immunol*. 2002;110(2 Suppl):S96-S106.

xxviii Rubber gloves: "born" - and now banished - at Johns Hopkins -- move addresses potentially fatal allergy to latex.

http://www.hopkinsmedicine.org/news/media/releases/rubber_gloves_born_and_now_banished_at_johns_hopkins.
_ Diakses 6 November 2015. d

xxix Edlich RF, Long WB, Gubler DK, et al. Dangers of cornstarch powder on medical gloves-seeking a solution. *Ann Plast Surg*. 2009; 63(1):822-826.

xxx UK Healthcare. University of Kentucky/UK HealthCare policy and procedure: latex-safe environment of care.

<http://www.hosp.uky.edu/policies/viewpolicy.asp?PolicyManual=10&PolicyID=966>. Diakses 6 November 2015.

xxxi Brodstone Memorial Hospital. Latex balloon guidelines.

<http://www.brodstonehospital.org/patients-and-visitors/latex-balloon-guidelines.html>. Diakses 6 November 2015.

xxxii Total Burn Care. Resident orientation manual: latex allergies.

http://www.totalburncare.com/orientation_latex_allergies.htm. Diakses 6 November 2015.

xxxiii Gibilisco PA. Clinical perils: latex allergy.

<http://www.icumed.com/media/84928/M1-1252%20Clinical%20Perils%20of%20Latex%20Rev.03PQ.pdf>. Diakses 6 November 2015.

xxxiv Brehler R, Kütting B. Natural rubber latex allergy- a problem of interdisciplinary concern in medicine.

<http://archinte.jamanetwork.com/article.aspx?articleid=647971>. Diakses 23 November 2014.

xxxv Brehler R, Kütting B. Natural rubber latex allergy- a problem of interdisciplinary concern in medicine.

<http://archinte.jamanetwork.com/article.aspx?articleid=647971>. Diakses 23 November 2014.

United States

Ansell Protective Products Inc.
111 Wood Avenue South
Suite 210
Iselin, NJ 08830, USA

EMEA

Ansell Healthcare Europe NV
Riverside Business Park
Blvd International, 55
1070 Brussels, Belgium

Australia

Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
Australia

Malaysia

Ansell Global Trading Center Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6,
63000 Cyberjaya, Malaysia

- xxxvi Surgical site infections. <https://www.vdh.virginia.gov/Epidemiology/Surveillance/HAI/ssi.htm>. Diakses 23 November 2015.
- xxxvii Kusachi S, Kashimura N, Konishi T, et al. Length of stay and cost for surgical site infection after abdominal and cardiac surgery in Japanese hospitals: multi-center surveillance. *Surg Infect (Larchmt)*. 2012;13(4):257-265.
- xxxviii Weber WP, Zwahlen M, Reck S, et al. Economic burden of surgical site infections at a European university hospital. *Infect Control Hosp Epidemiol*. 2008;29(7):623-629.
- xxxix Ranta PM, Ownby DR. A review of natural-rubber latex allergy in health care workers. *Clin Infect Dis*. 2004;38(2):252-256.
- xl Phillips VL, Goodrich MA, Sullivan TJ. Health care worker disability due to latex allergy and asthma: a cost analysis. *Am J Public Health*. 1999;89(7):1024-1028.
- xli Ranta PM, Ownby DR. A review of natural-rubber latex allergy in health care workers. *Clin Infect Dis*. 2004; 38(2):252–256.

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