

Ansell

**DOUBLE GLOVE
ONE DON. DONE**



GAMMEX® PI PLUS GLOVE-IN-GLOVE™ SYSTEM

Experience the convenience of the world's first pre-donned double gloving system. One don and you're done!

Features Ansell's softer and more durable next-generation polyisoprene (PI) gloves made with PI-KARE™ Technology that helps eliminate the use of standard chemical accelerators known to cause Type IV allergies and sensitivities.

Outer glove:


GAMMEX® Non-Latex PI Textured

- Textured fingers and palm for exceptional grip under dry and damp conditions
- DERMASHIELD™ Technology internal coating facilitates easy intra-operative changes of the top glove

Inner glove:

GAMMEX® Non-Latex PI Underglove

- Micro-thin underglove offers optimal comfort when double gloving
- Green color for easy glove breach identification
- DERMASHIELD™ Technology internal coating facilitates easy dry and damp hand donning



Reduces impact to the environment

Pre-donned outer and inner gloves, packed in one poly-pouch and one inner-wrap

Outer gloves come off easily

Top glove comes off easily for intra-operative glove changes

Pre-donned outer and inner gloves

Gloves are aligned, glove-to-glove, with minimal adjustments needed once donned



**Convenience of
two gloves-in-one**



**Saves half the
donning time**