

BIOCLEAN PERMAFLOW Sterile Cleanroom Pen S-BPFP

About

BioClean Permaflow Cleanroom Pens are sterile with fade and water resistant super permanent ink. With alcohol-resistant print on barrel for easy cleaning within a critical environment, they are ISO Class 4 & EU GMP Grade A compatible.

Specifications

MATERIAL:

Plastic barrel and lid, super permanent ink

LENGTH:

140mm (5.5")

COMPATIBILITY:

ISO Class 4 & EU GMP Grade A

CONSTRUCTION:

Plastic construction with super permanent ink.

COLOUR:

Black, Blue, Red or Green ink

STERILIZATION:

Gamma irradiation

Features

- Autoseal prevents drying
- Alcohol resistant print on barrel
- Processed to ensure ISO Class 4 & EU GMP Grade A compatibility
- Super permanent ink
- Quick drying
- Non-toxic


ALCOHOL RESISTANT PRINT ON BARREL

SUPER PERMANENT INK

FADE AND WATER RESISTANT

QUICK DRYING

NON-TOXIC



BioClean



Quality Standards

- Manufactured in a facility holding the ISO 9001:2015 certification
- Processed in an NEBB certified ISO Class 4 cleanroom environment



Sterilization

Method: Gamma irradiation
Minimum Dose: 25kGy
Sterility Assurance Level: 10^{-6}

Shelf Life & Storage

Five (5) years from date of manufacture.
Store in a dry, cool place (<40°C) away
from direct sunlight and fluorescent
light.

TO ORDER

Re-order code

Size

S-BPFP-BK1	Black
S-BPFP-BL1	Blue
S-BPFP-RD1*	Red
S-BPFP-GR1*	Green

Packing

One pen per sealed inner PE bag; 10 inner bags per sealed outer PE bag; 20 outer bags per lined carton (200 pens)

*Non-stock item, please contact for MOQ and lead times

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd

Ground Floor, 15 Kings
Court
Willie Snaith Road,
Newmarket,
Suffolk, CB8 7SG,
United Kingdom
Tel: +44 1638 663338

Ansell US

111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859

S-BPFP/1906/PDS22

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2020 Ansell Limited. All Rights Reserved.

