

# SANDEL® STAT-BLOC™ LINENS

O . R . T U R N O V E R S O L U T I O N S

## Disposable, antimicrobial covers to protect equipment and help provide a clean patient surface

### STAT-BLOC table sheets provide three layers of protection:

1. Quilted, breathable top layer allows for fluid absorption and minimises wrinkling
2. Strong concentration of Super Absorbent Polymers (SAP) for maximum fluid control
3. Impervious backing with proprietary antimicrobial technology protects the O.R. table from fluids

### STAT-BLOC™ Antimicrobial Arm Board Covers:

- Include dual straps for added safety
- Straps secured to prevent them from falling off and becoming contaminated

### STAT-BLOC™ Antimicrobial Headrest Covers:

- Fit over most standard head positioners
- Feature integrated flap for complete coverage

### STAT-BLOC™ Antimicrobial Stirrup Covers:

- Universal design to cover most stirrup styles


## KEY FEATURES & BENEFITS

- Antimicrobial linens are tested to be effective against MRSA<sup>1</sup>, E. coli<sup>1</sup>, CRE<sup>1</sup> and COVID-19<sup>2</sup>
- Antimicrobial linens are tested to show protection against bacterial migration to prevent colonisation of the linens<sup>3</sup>
- Helps reduce the amount of lint introduced to the operating room compared to standard reusable linens

## FEATURED TECHNOLOGIES



➔ [ansell.com/STAT-BLOC](https://ansell.com/STAT-BLOC)



## SANDEL®

### Anticipate the Unexpected

Developed in collaboration with clinical professionals, SANDEL products are designed to address key high-risk areas in the healthcare environment. These innovative safety solutions enable healthcare workers to do their job safely while improving patient outcomes.

# SANDEL® STAT-BLOC™ LINENS

O . R . T U R N O V E R S O L U T I O N S

## Disposable draw sheets to aid patient transfer and help ensure cleanliness

### SANDEL Draw Sheets:

- Hydrophilic to allow fluid to easily pass through and be absorbed by the table sheet
- 450 kg. weight limit in wet and dry conditions<sup>4</sup>

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
AS-30-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Table Sheet (76 x 101 cm)	50/CS	50
AS-60-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Table Sheet (152 x 101 cm)	50/CS	50
AS-92-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Table Sheet (228 x 101 cm)	20/CS	20
AS-92-QZF	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Table Sheet (228 x 101 cm); Fitted	20/CS	20
AS-100-QZ	STAT-BLOC™ Table Sheet	Disposable, Antimicrobial, Absorbent Table Sheet (254 x 101 cm)	20/CS	20
ABC-10	STAT-BLOC™ Arm Board Cover	Disposable, Antimicrobial Arm Board Cover (33 x 71 cm) with 2 Arm Straps	20/Bag; 10 Bags/CS	200
HR-10	STAT-BLOC™ Headrest Cover	Disposable, Antimicrobial Headrest Cover (38 x 38 cm)	25/Bag; 10 Bags/CS	250
SC-01	STAT-BLOC™ Stirrup Cover	Disposable, Antimicrobial Stirrup Cover (38 x 38 cm)	20/Bag; 10 Bags/CS	200
LFS-66	Draw Sheet	Disposable Draw Sheet (86 x 168 cm)	5/Bag; 10 Bags/CS	50



AS-92-AZ



HR-10



ABC-10



SC-01



LFS-66

### REFERENCES

- 1 Data on file: Antimicrobial Protection Using STAT-BLOC™ Table Sheets, Environmental Protection for Patient Safety
- 2 Data on file: Summary of Testing on Virucidal Properties of STAT-BLOC™ Table Sheets when Challenged with Coronavirus
- 3 Data on file: Clinical Relevance Testing Report: An evaluation of bacterial migration from a contaminated surface to a patient surrogate with the use of a variety of table covers
- 4 Data on file: Testing Report: Lift Strength of SANDEL® Lift Sheet

### ✦ Contact your Ansell representative to order or for more information.

**North America**  
 US Tel: 866-764-3327  
 818-534-2500  
 Email: sandel-cs@ansell.com

**Latin America & Caribbean**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
 Brazil Tel: +55-11-3356-3100  
 Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0)2 528 74 00  
 Email: info@ansell.eu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com  
**Japan**  
 Tel: +81 (0)3 5549 8151  
 Email: info.medical.jp@ansell.com

**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com