

# SANDEL® STAT-BLOC™ ANTIMICROBIAL LINENS

O . R . T U R N O V E R S O L U T I O N S

## Disposable, antimicrobial table covers, arm board covers and headrest covers to help prevent cross-contamination

### STAT-BLOC table sheets provide three layers of protection:

1. Quilted, breathable top layer allows for fluid absorption and minimizes wrinkling to help prevent skin breakdown
2. Strong concentration of Super Absorbent Polymers (SAP) for maximum fluid control
3. Impervious backing with proprietary antimicrobial technology offers complete protection against bacterial migration<sup>1</sup>

STAT-BLOC arm board covers feature dual, secured straps for added safety and to prevent straps from falling to the floor and becoming contaminated. STAT-BLOC headrest covers feature an integrated flap for full coverage.

### KEY FEATURES & BENEFITS


- AMT™ Antimicrobial Technology is proven 99.9% effective against MRSA, E. coli and CRE and 99% effective against COVID-19<sup>1</sup>
- Shown in testing to be more absorbent than leading table sheets and absorb 3x faster than the leading table sheet<sup>1</sup>
- Table sheets are made of 49% renewable bio-based material which is carbon neutral when incinerated<sup>1</sup>
- Orange color is a visible indicator of the antimicrobial and provides separation from blue sterile drapes
- Less lint than standard reusable linens to help prevent cross-contamination
- Lift sheet and armboard covers come folded inside the table sheet

### FEATURED TECHNOLOGIES



[ansell.com/STAT-BLOC](https://ansell.com/STAT-BLOC)

<sup>1</sup> Data on file.



## SANDEL® Anticipate the Unexpected

Through collaboration with clinical professionals, SANDEL products are designed to address key high-risk issues affecting healthcare workers. Its signature high-visibility orange products and premium antimicrobial technology equip caregivers in the healthcare environment to do their job safely while improving outcomes.

# SANDEL® STAT-BLOC™ ANTIMICROBIAL LINENS

O . R . T U R N O V E R S O L U T I O N S

ITEM #	PRODUCT NAME	DESCRIPTION	CASE QTY.
SP-02	STAT-BLOC Linen Set (ORANGE)	Antimicrobial, Absorbent Table Sheet, Lift Sheet, Arm Board Covers	20
SP-03	STAT-BLOC Linen Set (ORANGE)	Antimicrobial, Absorbent Table Sheet, Lift Sheet, Arm Board Covers, Headrest Cover	20
SP-05	STAT-BLOC Linen Set (ORANGE) with Strap	Antimicrobial, Absorbent Table Sheet, Lift Sheet, Arm Board Covers and Strap	15
SP-06	STAT-BLOC Linen Set (ORANGE) with Strap and Microfiber Mop	Antimicrobial, Absorbent Table Sheet, Lift Sheet, Arm Board Covers, Strap, and Microfiber Mop	15

Linen sets can be included in custom turnover packs with other components such as bags, patient positioning straps, mops and wipes.



STAT-BLOC Table Sheet (Folded)



STAT-BLOC Headrest Cover



STAT-BLOC Armboard Cover

➤ **Contact your Ansell representative for ordering or more information.**

**North America**  
 US Tel: 866-764-3327  
 818-534-2500  
 Email: sandel-cs@ansell.com

**Latin America & Caribbean**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
 Brazil Tel: +55-11-3356-3100  
 Brazil Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0)2 528 74 00  
 Email: info@ansellLeu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com