

## Supreno® Plus 93-753

Durable Nitrile Exam Glove with Extended Cuff and Advanced Barrier Protection

# Nitrile disposable gloves offering extended chemical splash protection

- Strong formulation: MICROFLEX<sup>®</sup> Supreno<sup>®</sup> Plus 93-753 gloves' robust nitrile formulation makes them more durable and resistant to rips, snags and tears
- **Real-world performance:** These robust disposable gloves boast EN ISO Type B chemical splash protection, making them suitable for professionals in high-risk environments
- Low AQL: MICROFLEX<sup>®</sup> Supreno<sup>®</sup> Plus 93-753 gloves' low 0.65 AQL reinforces their barrier integrity
- Extended protection: Featuring an extended cuff, these gloves offer forearm and wrist protection for workers



#### Industries

- Automotive
- Automotive Aftermarket
- Life Sciences

#### **Recommended For**

- Inspection, selecting, checking parts
- Assembly and inspection of components
- Picking and fastening components
- Equipment repair and maintenance
- General Purpose Auto Aftermarket
- Tire rotation & changes
- Animal care feeding and cleanup
- Maintaining lab furniture and equipment
- Sample taking and processing
- Administering drugs
- Raw material sample collection
- Extra protection over the wrist and arm



### **TECHNICAL DATA SHEET**

	Product Information
External Glove Surface	Textured Fingers
Audit Standards	ISO 13485, ISO 9001
Packaging Overview	50 gloves per dispenser; 10 dispensers per case; 500 gloves per case;
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country Of Origin	Malaysia
Product Segmentation	High Risk
Antistatic	Yes

#### **PHYSICAL PROPERTIES**

Typical Values		Testing Method
295 / 11.6		ASTM D3767,EN 420
0.65 AQL		ASTM D5151,EN 374-2
0.14 / 5.5		ASTM D3767,EN 420
0.22 / 8.7		ASTM D3767,EN 420
BEFORE AGING	AFTER AGING	
32	27	ASTM D412 & D573
550	500	ASTM D412
15	12	EN 455-2
	295 / 0.65 0.14 , 0.22 , <b>BEFORE AGING</b> 32 550	295 / 11.6 0.65 AQL 0.14 / 5.5 0.22 / 8.7 BEFORE AGING AFTER AGING 32 27 550 500

#### **ORDERING INFORMATION**

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	93753070	93753080	93753090	93753100

#### **Performance Standards and Regulatory Compliance**



#### For additional information visit us at www.ansell.com, or call us at

latam.medical@ansell.com MX: +52 442 296 2050 latam.medical@ansell. CO: +57 1 288 3247 latam.medical@ansell.

#### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV T: +32 (0) 2 528 74 00 F: +32 (0) 2 528 74 01

#### Asia Pacific Region MY +603 8310 6688 apac.medical@ansell.com CN: +86 21 38275000

infochir JP: +813 5549 8151 al.jp@ansell.com info.m

North America Region Ansell Healthcare Products LLC T: +1 800 800 0444 F: +1 800 800 0445 Latin America & Caribbean Region BR +55 (11) 3356-3100 Australia Ansell Limited T: +61 1800 337 041 F: +61 1800 803 578 Russia

Ansell PYC Ten. +7 495 258 13 16

Ansell, <sup>®</sup> and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

ell.com

ell.com

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

