

# GAMMEX<sup>®</sup> PI Ergo-Fit<sup>™</sup>

S U R G I C A L   G L O V E S

## Polyisoprene surgical gloves, with back-of-hand and finger relief zones for added comfort

- **Elevated comfort:** Thanks to their ridged back-of-hand and finger relief zones, GAMMEX<sup>®</sup> PI Ergo-Fit sterile surgical gloves feel more spacious and require 20% less muscle effort\*, resulting in an optimized range of motion and barehand-like comfort
- **Reduced allergy risks:** These polyisoprene non-latex disposable gloves are free from diphenylguanidine (DPG) and cetylpyridium chloride (CPC), eliminating the risk of latex Type I allergies and minimizing chemical Type IV allergies and sensitivities
- **Facilitated donning:** GAMMEX<sup>®</sup> PI Ergo-Fit disposable surgical gloves features DERMASHIELD<sup>™</sup> Technology, an internal coating that allows for fast and easy donning in dry and damp conditions  
\*compared to standard PI surgical gloves lacking this feature



### Recommended For

- All surgeries, especially for long, sustained procedures
- Use as an outer glove when double gloving

### FEATURED TECHNOLOGIES



### GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> PI Ergo-Fit<sup>™</sup>

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White
Thickness	Standard
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD <sup>™</sup> Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, EN 374-1, 2, 3, EN 420:2003 + A1:2009, EN 455 1-4, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340064055	340064060	340064065	340064070	340064075	340064080	340064085	340064090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com