

Validated sterile nitrile RABS/isolator sleeve featuring low particle levels and superior ESD properties

- **Assured quality:** The BioClean™ Sterile Nitrile RABS/Isolator Sleeve GSL36NITPP26 is 100% visually inspected and triple bagged to ensure high-quality arm protection
- **Enhanced versatility:** As this isolator sleeve is compatible with autoclaving, VHP and IPA, it is suitable for three distinct sterilization methods
- **Increased practicality:** The BioClean™ Sterile Nitrile RABS/Isolator Sleeve GSL36NITPP26 is also compatible with a broad range of cuff systems
- **Added features:** In addition, this sleeve has extremely low particulate levels and excellent ESD properties



KEY FEATURES & BENEFITS

- **Triple bagged and 100% visually inspected:** Excellent arm protection
- **Versatile RABS/Isolator sleeves:** Compatible with a variety of cuff systems
- **Autoclave compatibility:** Suitable for pressurized, high-temperature conditions

Industries

- Pharmaceutical Manufacturing
- Laboratory Research and Development
- Raw Material and Preparation
- Biotechnology Manufacturing
- Sterile Processing Department
- Semiconductor manufacturing





GSL36NITPP26

Validated Sterile Nitrile RABS / Isolator Sleeve

TECHNICAL DATA SHEET

Product Reference	GSL36NITPP26
Material	Nitrile
Color	White
External Glove Surface	Smooth
Cuff Style	Beaded
Audit Standards	ISO 14001, ISO 9001:2008, NEBB Certified Cleanrooms
Packaging Overview	One piece per sealed inner PE bag; one inner bag per sealed second inner PE bag; one second inner bag per sealed outer PE bag; 20 outer bags per lined inner white Correx polyethylene box (20 pieces).
Storage Instructions	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country Of Origin	Malaysia
Temperature Range	From -10°C to 120°C
Sterilization Method	GAMMA irradiation (25 kGy)
Cleanroom Class	Class 10/ISO 4 & EU GMP Grade A
Shelf Life	Three (3) years from date of manufacture.
Protein Level	N/A: contains no natural rubber latex
Antistatic	Yes
Standards	Category III, CE 0493, EC 1935/2004, ISO 11137-1:2006, EN 388:2016 +A1:2018, EN 420:2003 + A1:2009, EN 455 Part 2, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 421:2010



GSL36NITPP26

Validated Sterile Nitrile RABS / Isolator Sleeve

PHYSICAL PROPERTIES

	Typical Values	Testing Method
Typical Length (mm/in)	660/26	EN ISO 21420
Port Size	12"	
Target Single Wall Cuff Thickness (mm/mil)	0.50/20	EN 455-2

Technology



NO RUBBER LATEX

Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890



EN ISO 374-1:2016

EN ISO 374-5



VIRUS

EN 421



AJKNOTS

EN588: 2016



4102X

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

