

Previously known as: MICROGARD® 2500 Apron - Model 213

Durable microporous laminate protective apron, keeping wearers safe against biological agents.

- **Specialized protection:** This AlphaTec® protective apron achieves the highest classifications of protection against biological agents, in accordance with EN 14126 for penetration of blood, body fluids and blood-borne pathogens
- **Improved durability:** Its unique microporous polyethylene laminate nonwoven material offers excellent strength and durability, for personal protection with greater longevity
- **Reduced contamination risks:** As this safety apron's material is also low-linting by design, the chances of contamination in critical areas are minimized



Industries

- Life Sciences
- Industrial
- Chemical
- Law Enforcement

Applications


- Transferring liquids and solids
- Unexpected leakages spills or other releases
- Plant disinfection and sanitization
- Emergency Response
- Spill or leakage cleanup
- Maintenance
- Contaminated Environments

KEY FEATURES & BENEFITS


- **Protective barrier:** High class of protection against biological agents*
- **Unique microporous laminate formulation:** Enhanced strength and durability
- **Low-linting material:** Reduced contamination risk in critical areas

*in accordance with EN 14126, for penetration of blood, body fluids and blood-borne pathogens

Material Diagram



Seam Diagram



Performance Standards



EN ISO 13688:2013



Product Information

Product Information	
Available Sizes	One Size
Color	White
Country Of Origin	China
Packaging Overview	50 every compressed bag
Product Material	Microporous polyethylene Laminate non-woven
Product Reference	2500 Apron - Model 213
Seam Type	Stitched
Shelf Life	5 years
Style Number	W25S213

Suit components & accessories

Suit Components and Accessories

Design Features

Apron with tie fastening to waist and 100cm long tie fastening

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00

North America Region

Ansell Healthcare Products LLC
T: +1 800 800 0444

Australia

Ansell Limited
T: +61 1800 337 041

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

