

## Sterile disposable latex gloves designed for cleanroom use, offering both flexibility and comfort

- **Heightened comfort and flexibility:** BioClean™ Alpha AL300 protective gloves are anatomically designed and made from sterile natural rubber latex, ensuring both comfort and flexibility
- **Improved durability:** These disposable powder-free latex gloves boast a reinforced cuff design, boosting strength and stability while reducing instances of glove tearing when donning
- **Enhanced usability:** They are also designed with EasyOn technology, enabling double gloving with ease
- **Assured cleanliness:** BioClean™ Alpha AL300 gloves come in non-particulating EasyTear packaging, for genuinely cleanroom-compatible hand protection



### KEY FEATURES & BENEFITS

- **Hand-specific sterile latex design:** Superb comfort and flexibility
- **Reinforced cuff:** Fewer instances of tearing when donning
- **Smooth grip:** For more confident, secure handling

### Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D
- Sterile Processing Department



# Alpha AL300

Sterile Disposable Latex Cleanroom Glove

## TECHNICAL DATA SHEET

Product Information	
Material	Natural Rubber Latex
Color	Natural
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	CE 0493, EN ISO 21420:2020, EN 421:2010, EN 455 Part 2, EN ISO 374-1:2016, EN ISO 374-5:2016, Food Contact, Category III, UKCA
Packaging	One pair per inner PE wallet; one wrap per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs)  <b>More sustainable packaging:</b> These sterile latex gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Malaysia
Available sizes	6, 6.5, 7, 7.5, 8, 8.5, 9
External Glove Surface	Textured Fingers
Internal Glove Surface	Chlorinated
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 <sup>-6</sup>
Cleanroom Class	Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	Five (5) years from date of manufacture.
Protein Level	50 µg/g or less of total extractable protein
Anti-static	No



# Alpha AL300

Sterile Disposable Latex Cleanroom Glove

Physical Properties								Testing Method
Sizes	6	6.5	7	7.5	8	8.5	9	
Typical Length (mm/in)	295 / 11.61							EN ISO 21420
Palm Width (mm/in)	77/3	83/3.3	89/3.5	95/3.7	102/4	109/4.3	114/4.5	
Freedom from Holes	1.5 AQL Performance Level 2							EN 374-2
Typical Particle Count $\geq 0.5\mu\text{m}$ (counts / $\text{cm}^2$ )	<3500							IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.20 / 7.87							EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.24 / 9.44							EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.18 / 7.08							EN 455-2
Ultimate tensile strength (MPa) During Aging	Min. 14							ASTM D412-06a
Force at Break (N) During Aging	$\geq 9\text{N}$							EN 455-2

## IONIC CONTENT

Concentration in $\mu\text{g}/\text{cm}^2$	Typical	Concentration in $\mu\text{g}/\text{cm}^2$	Typical
Ammonium	< 0.05	Nitrate	< 0.50
Bromide	Available on File	Nitrite	< 0.01
Calcium	< 1.00	Phosphate	< 0.01
Chloride	< 2.50	Potassium	< 0.08
Fluoride	< 0.01	Sodium	< 0.08
Lithium	< 0.01	Sulphate	< 0.20
Magnesium	< 0.05	Zinc	< 1.20

## ORDERING INFORMATION

SIZE	6	6.5	7	7.5	8	8.5	9
REORDER NO.	AL300-60	AL300-65	AL300-70	AL300-75	AL300-80	AL300-85	AL300-90

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

