

# GAMMEX® PI Plus Glove-in-Glove™ System

S U R G I C A L   G L O V E S

## Polyisoprene inner and outer surgical gloves, pre-donned and ready for double gloving

- **Effortless double gloving:** The GAMMEX® PI Plus Glove-in-Glove™ System is a unique two-in-one system with both inner and outer polyisoprene surgical gloves already pre-donned and aligned at the fingertips, enabling wearers to double glove in one don and save half the donning time
- **Sustainable packaging:** Being pre-donned, the inner and outer sterile gloves are packed in one inner wrap for a 50% reduction of paper
- **Advanced allergy prevention:** These non-latex gloves present no risk of latex allergies and do not contain chemical accelerators known to cause allergic contact dermatitis
- **Ergonomic comfort:** Certified by United States Ergonomics for superb wearer comfort and performance
- **Effortless donning and doffing:** Ansell's donning aid technology facilitates easy donning and doffing of the outer glove during intraoperative changes



## Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach

## FEATURED TECHNOLOGIES



## GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System

SURGICAL GLOVES

## PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	White, Green
Thickness	Standard
Cuff Style	Beaded with SUREFIT <sup>™</sup> Technology
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated with DERMASHIELD <sup>™</sup> Technology
Allergy Prevention	Latex (Type I)
Grip Level	High
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	2 pairs per pouch; 50 pairs per box (25 pouches); 4 boxes per carton/case; 200 pairs per carton/case (100 pouches)
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Xanthates, Zinc Dialkyldithiocarbamate

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340082055	340082060	340082065	340082070	340082075	340082080	340082085	340082090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com