

Previously known as: MICROGARD® 2500 Apron - Model 213

### Durable microporous laminate protective apron, keeping wearers safe against biological agents.

- **Specialized protection:** This AlphaTec® protective apron achieves the highest classifications of protection against biological agents, in accordance with EN 14126 for penetration of blood, body fluids and blood-borne pathogens
- **Improved durability:** Its unique microporous polyethylene laminate nonwoven material offers excellent strength and durability, for personal protection with greater longevity
- **Reduced contamination risks:** As this safety apron's material is also low-linting by design, the chances of contamination in critical areas are minimized



## Industries

- Life Sciences
- Industrial
- Chemical
- Law Enforcement

## Applications


- Transferring liquids and solids
- Unexpected leakages spills or other releases
- Plant disinfection and sanitization
- Emergency Response
- Spill or leakage cleanup
- Maintenance
- Contaminated Environments

## KEY FEATURES & BENEFITS


- **Protective barrier:** High class of protection against biological agents\*
- **Unique microporous laminate formulation:** Enhanced strength and durability
- **Low-linting material:** Reduced contamination risk in critical areas

\*in accordance with EN 14126, for penetration of blood, body fluids and blood-borne pathogens

## Material Diagram



## Seam Diagram



## Performance Standards



EN ISO 13688:2013



## Product Information

Product Information	
Available Sizes	One Size
Color	White
Country Of Origin	China
Packaging Overview	50 every compressed bag
Product Material	Microporous polyethylene Laminate non-woven
Product Reference	2500 Apron - Model 213
Seam Type	Stitched
Shelf Life	5 years
Style Number	W25S213

## Suit components & accessories

### Suit Components and Accessories

#### Design Features

Apron with tie fastening to waist and 100cm long tie fastening

**For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at**

**Europe, Middle East & Africa Region**

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00

**North America Region**

Ansell Healthcare Products LLC  
T: +1 800 800 0444

**Australia**

Ansell Limited  
T: +61 1800 337 041

**Asia Pacific Region**

Ansell Global Trading Center  
T: +603 8310 6688

**Latin America & Caribbean Region**

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

