




Ultimate BUPS

Sterile Disposable Neoprene
(Polychloroprene) Cleanroom Glove

Sterile polychloroprene cleanroom gloves, offering users chemical resistance and tactility

- **Specialized protection:** BioClean™ Ultimate BUPS cleanroom gloves are tested against ASTM D6978 and EN standards, guarding against chemotherapy drugs and various chemicals
- **Assured component integrity:** These latex-free gloves' electrostatic-dissipative (ESD) properties ensure greater component protection against electrostatic discharge
- **Reduced contamination risks:** Their non-particulating EasyTear™ packaging is designed to minimize contamination risks
- **Enhanced features:** These powder-free cleanroom gloves are manufactured with a beaded cuff, enhancing their strength and facilitating double-donning
- **Heightened tactility and comfort:** BioClean™ Ultimate BUPS polychloroprene gloves are also comfortable, with good tactility levels for engaging in precision work



KEY FEATURES & BENEFITS

- **ASTM and EN standards:** Chemical and chemotherapy drug resistance
- **ESD properties:** Component integrity protection
- **Non-particulating EasyTear packaging:** Fewer contamination risks

Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D
- Sterile Processing Department

WWW.ANSELL.COM





Ultimate BUPS

Sterile Disposable Neoprene
(Polychloroprene) Cleanroom Glove

TECHNICAL DATA SHEET

| Product Information | |
|--|---|
| Material | Neoprene (Polychloroprene) |
| Color | Natural |
| Shape | Hand Specific |
| Cuff | Beaded |
| Manufacturing/QMS Audit Standards | ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms |
| Regulatory/Standard Compliance | ASTM D6978, CE 0598, EN ISO 21420:2020, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Category III |
| Packaging | One pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs) |
| Storage | Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light. |
| Country of Origin | Indonesia |
| Available sizes | 6, 6.5, 7, 7.5, 8, 8.5, 9 |
| Powder Content | Powder-Free |
| External Glove Surface | Textured Fingers and Palm |
| Internal Glove Surface | Polymer Coated |
| Sterilization Method | GAMMA irradiation (25 kGy) |
| Sterilization Minimum Dose | 25kGy |
| Sterility Assurance Level | 10 ⁻⁶ |
| Cleanroom Class | Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms |
| Shelf Life | Five (5) years from date of manufacture. |
| Tested for use with Chemotherapy Drugs | Yes |
| Protein Level | N/A: contains no natural rubber latex |
| Anti-static | Yes |



Ultimate BUPS

Sterile Disposable Neoprene
(Polychloroprene) Cleanroom Glove

| Physical Properties | | | | | | | | Testing Method |
|--|------------------------------|--------|--------|--------|-------|---------|---------|-----------------|
| Sizes | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 | |
| Typical Length (mm/in) | 300 / 12 | | | | | | | EN 420 |
| Palm Width (mm/in) | 77/3 | 83/3.2 | 89/3.5 | 95/3.7 | 102/4 | 108/4.2 | 115/4.5 | |
| Freedom from Holes | 0.65 AQL Performance Level 3 | | | | | | | EN 374-2 |
| Typical Particle Count ≥0.5µm (counts / cm ²) | <1300 | | | | | | | IEST-RP-CC005.4 |
| Target Single Wall Palm Thickness (mm/mil) | 0.11 / 4.33 | | | | | | | EN 455-2 |
| Target Single Wall Finger Thickness (mm/mil) | 0.14 / 5.51 | | | | | | | EN 455-2 |
| Target Single Wall Cuff Thickness (mm/mil) | 0.09 / 3.54 | | | | | | | EN 455-2 |
| Force at Break (N) Before Aging | ≥ 9 N | | | | | | | EN 455-2 |
| Force at Break (N) During Aging | ≥ 9 N | | | | | | | EN 455-2 |

IONIC CONTENT

| Concentration in µg/cm ² | Typical | Concentration in µg/cm ² | Typical |
|-------------------------------------|--------------|-------------------------------------|--------------|
| Ammonium | 0.005 | Nitrate | 0.795 |
| Bromide | Not Detected | Nitrite | Not Detected |
| Calcium | 0.582 | Phosphate | Not Detected |
| Chloride | 0.494 | Potassium | 0.378 |
| Fluoride | Not Detected | Sodium | 0.123 |
| Lithium | Not Detected | Sulphate | 0.041 |
| Magnesium | 0.025 | Zinc | Not Tested |

ORDERING INFORMATION

| SIZE | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 |
|-------------|---------|---------|---------|---------|---------|---------|---------|
| REORDER NO. | BUPS-60 | BUPS-65 | BUPS-70 | BUPS-75 | BUPS-80 | BUPS-85 | BUPS-90 |

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

