

# GAMMEX® Non-Latex PI Underglove

S U R G I C A L   G L O V E S

## Soft, skin-friendly polyisoprene surgical gloves, for heightened tactile sensitivity

- **Superior comfort and protection:** Ansell's proprietary soft yet strong formulation delivers heightened comfort and hand protection when double gloving
- **Easy breach detection:** Green in color, these sterile gloves enable fast and easy breach detection when worn with a cream or white outer glove
- **Optimal allergy protection:** Being non-latex, these polyisoprene surgical gloves eliminate the risk of developing a latex allergy while their skin-friendly formulation helps minimize sensitization to developing a Type IV chemical allergy as chemical accelerators known to cause skin allergies have been eliminated
- **Effortless donning and doffing:** A unique internal coating enables easy donning and doffing under dry or damp hand conditions



## Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures

## FEATURED TECHNOLOGIES



## GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Non-Latex PI Underglove

SURGICAL GLOVES

## PRODUCT INFORMATION

<b>Material</b>	PI (Polyisoprene)
<b>Powder Content</b>	Powder-Free
<b>Color</b>	Green
<b>Thickness</b>	Micro
<b>Cuff Style</b>	Beaded with SUREFIT <sup>™</sup> Technology
<b>External Glove Surface</b>	Smooth with Micro-Textured Finish
<b>Internal Glove Surface</b>	Polymer Coated with DERMASHIELD <sup>™</sup> Technology
<b>Allergy Prevention</b>	Latex (Type I)
<b>Grip Level</b>	Low
<b>Double Gloving Recommendation</b>	Underglove
<b>Tested For Use With Chemotherapy Drugs</b>	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
<b>Freedom from Holes (Inspection level I)</b>	0.65 AQL
<b>Sterile</b>	Yes
<b>Product Standards</b>	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
<b>Quality / Environmental Standards</b>	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
<b>Packaging</b>	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
<b>Shelf Life</b>	3 years
<b>Storage Instructions</b>	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
<b>Vulcanization Chemical Accelerators</b>	Xanthates, Zinc Dialkyldithiocarbamate

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
<b>Product Code</b>	340058055	340058060	340058065	340058070	340058075	340058080	340058085	340058090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com