

# GAMMEX® PI Plus Glove-in-Glove™ System

S U R G I C A L   G L O V E S

## Polyisoprene inner and outer surgical gloves, pre-donned and ready for double gloving

- **Effortless double gloving:** The GAMMEX® PI Plus Glove-in-Glove™ System is a unique two-in-one system with both inner and outer polyisoprene surgical gloves already pre-donned and aligned at the fingertips, enabling wearers to double glove in one don and save half the donning time
- **Sustainable packaging:** Being pre-donned, the inner and outer sterile gloves are packed in one inner wrap for a 50% reduction of paper
- **Advanced allergy prevention:** These non-latex gloves present no risk of latex allergies and do not contain chemical accelerators known to cause allergic contact dermatitis
- **Ergonomic comfort:** Certified by United States Ergonomics for superb wearer comfort and performance
- **Effortless donning and doffing:** Ansell's donning aid technology facilitates easy donning and doffing of the outer glove during intraoperative changes



### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Easy identification of a glove breach

### FEATURED TECHNOLOGIES



## GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> PI Plus Glove-in-Glove<sup>™</sup> System

SURGICAL GLOVES

## PRODUCT INFORMATION

|   |   |
|---|---|
| Material                                | PI (Polyisoprene)   |
| Powder Content                          | Powder-Free   |
| Color                                   | White, Green  |
| Thickness                               | Standard  |
| Cuff Style                              | Beaded with SUREFIT <sup>™</sup> Technology   |
| External Glove Surface                  | Textured Fingers and Palm   |
| Internal Glove Surface                  | Polymer Coated with DERMASHIELD <sup>™</sup> Technology   |
| Allergy Prevention                      | Latex (Type I)  |
| Grip Level                              | High  |
| Tested For Use With Chemotherapy Drugs  | Yes, in accordance with ASTM D6978 and US FDA cleared   |
| Freedom from Holes (Inspection level I) | 0.65 AQL  |
| Sterile                                 | Yes   |
| Product Standards                       | AS/NZS 4179, ASTM D3577, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107 |
| Quality / Environmental Standards       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001  |
| Packaging                               | 2 pairs per pouch; 50 pairs per box (25 pouches); 4 boxes per carton/case; 200 pairs per carton/case (100 pouches)                                    |
| Shelf Life                              | 3 years   |
| Storage Instructions                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.  |
| Vulcanization Chemical Accelerators     | Xanthates, Zinc Dialkyldithiocarbamate  |

## ORDERING INFORMATION

| Size         | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Product Code | 340082055 | 340082060 | 340082065 | 340082070 | 340082075 | 340082080 | 340082085 | 340082090 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com