

## Ambidextrous non-sterile latex cleanroom gloves, for added convenience and comfort

- **Heightened flexibility and comfort:** Thanks to their soft, stretchy film, BioClean™ Legacy BLA2 non-sterile PPE gloves combine excellent comfort levels with heightened flexibility
- **Optimized cleanliness:** Their low particle count and minimal ionic content make them ideal for cleanroom applications
- **Increased practicality:** Ambidextrous by design, BioClean™ Legacy BLA2 safety gloves are convenient to wear on either hand, while their smooth surface enables effortless double-donning

*CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.*



### KEY FEATURES & BENEFITS

- **Soft, stretchy film:** For safety gloves offering comfort and flexibility
- **Textured surface and beaded cuff:** Enhanced grip and ease of donning
- **Clean processed:** Low particle count and ionic content\*

\*suitable for cleanroom applications

### Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D

## TECHNICAL DATA SHEET

Product Information	
Material	Natural Rubber Latex
Color	Natural
Shape	Ambidextrous
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	EN 421:2010, EN 455 Part 2, EN ISO 374-1:2016, EN ISO 374-5:2016, Food Contact, Category III, UKCA
Packaging	100 pieces per sealed inner PE bag; one inner PE bag per sealed outer PE bag; 10 outer bags per lined carton (1000 pieces). <b>More sustainable packaging:</b> These non-sterile latex gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases. *Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.
Storage	Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.
Country of Origin	Malaysia
Available sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10), 2XL (10.5-11)
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Chlorinated
Cleanroom Class	Class 10/ISO 4
Shelf Life	Five (5) years from date of manufacture.
Tested for use with Chemotherapy Drugs	No
Protein Level	50 µg/g or less of total extractable protein
Anti-static	No



# Legacy BLA2

## Non-sterile Disposable Latex Cleanroom Glove

Physical Properties							Testing Method
Sizes	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	2XL (10.5-11)	EN ISO 21420
Typical Length (mm/in)	300 / 12						
Palm Width (mm/in)	77/3	83/3.3	95/3.7	106/4.2	112/4.4	118/4.6	
Freedom from Holes	1.5 AQL Performance Level 2						EN 374-2
Typical Particle Count $\geq 0.5\mu\text{m}$ (counts / $\text{cm}^2$ )	<1500						IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.17 / 6.69						EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.20 / 7.87						EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.11 / 4.33						EN 455-2
Force at Break (N) During Aging	$\geq 6$ N						EN 455-2

## IONIC CONTENT

Concentration in $\mu\text{g}/\text{cm}^2$	Typical	Concentration in $\mu\text{g}/\text{cm}^2$	Typical
Ammonium	Not Tested	Nitrate	0.67
Bromide	Not Tested	Nitrite	Not Tested
Calcium	0.37	Phosphate	0.07
Chloride	0.97	Potassium	0.03
Fluoride	Not Detected	Sodium	0.02
Lithium	Not Tested	Sulphate	Not Detected
Magnesium	0.003	Zinc	Not Detected

## ORDERING INFORMATION

SIZE	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	2XL (10.5-11)
REORDER NO.	BLA200-XS	BLA201-S	BLA202-M	BLA203-L	BLA204-XL	BLA205-XXL

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

