

Non-sterile, tie-on face mask with advanced filtration efficiency (single-packed)

- **Specialized protection:** This tie-on face mask offers high particle and bacterial filtration and is tested against ASTM F2101 standards for Bacterial Filtration Efficiency (BFE)
- **Enhanced features and comfort:** The BioClean™ Non Sterile Tie on Facemask MTA210-2 is equipped with a fully enclosed, malleable noseband, which enhances fit and comfort, while its ultrasonic seal leads to a stronger bond around the tie-tapes and enables secure fastening
- **Reduced contamination risks:** This non-sterile face mask is processed in a NEBB-certified ISO Class 4 cleanroom environment, minimizing potential contamination risks

Key Features and Benefits

- **Contamination control:** High-level Bacterial Filtration Efficiency (BFE)
- **Cleanroom-compatible materials:** Reduced contamination risks
- **Ultrasonically sealed edges:** Tie-tapes with a reliable, strong bond

Industries

- Production and Manufacturing
- Controlled and Critical Environments
- Laboratory and Research





BioClean™ Non Sterile Tie on Facemask MTA210-2

TECHNICAL DATA SHEET

PRODUCT INFORMATION

| | |
|-----------------------------|--|
| Material | Polypropylene Polyethylene |
| Audit Standards | Manufacturing QMS Audit Standards ISO 9001 |
| Standards | ASTM F2101, ISO 2859, ISO 9001:2015 |
| Packaging Overview | One piece per sealed inner PE bag; 50 inner bags per sealed outer PE bag; 10 outer bags per lined carton (500 pieces). |
| Country Of Origin | China |
| Storage Instructions | Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light. |
| Shelf Life | Five (5) years from date of manufacture. |

PHYSICAL PROPERTIES

| | |
|---------------------------|--|
| Color | White |
| String Type | Ties |
| Overall Dimensions | WIDTH (± 5mm) 210mm LENGTH (± 5mm) 95mm |
| Inner Facing Layer | Non-woven polypropylene/polyethylene (hygroscopic) |
| Filtration Layer | Meltblown polypropylene |
| Outer Layer | Non-woven polypropylene polyethylene hydrophobic |
| Fastenings | Non-woven polypropylene ties |
| Noseband | Plastic coated steel |

PERFORMANCE INFORMATION

| | |
|--|--|
| Cleanroom Class | Class 10/ISO 4 |
| Particle Filtration Efficiency | 99.9% @0.1 micron |
| Bacterial Filtration Efficiency | >99.9% @0.3 micron |
| Viral Filtration Efficiency | 99.8% @ 3.0 micron |
| Differential Pressure | 6.9mm H ₂ O/cm ² |

ORDERING INFORMATION

| | | |
|----------|-------------|---------------------------|
| | SIZE | Individual Packing |
| MTA210-2 | REORDER NO. | MTA210-2 |

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

