




# Ultimate BUPS

Sterile Disposable Neoprene  
(Polychloroprene) Cleanroom Glove

## Sterile polychloroprene cleanroom gloves, offering users chemical resistance and tactility

- **Specialized protection:** BioClean™ Ultimate BUPS cleanroom gloves are tested against ASTM D6978 and EN standards, guarding against chemotherapy drugs and various chemicals
- **Assured component integrity:** These latex-free gloves' electrostatic-dissipative (ESD) properties ensure greater component protection against electrostatic discharge
- **Reduced contamination risks:** Their non-particulating EasyTear™ packaging is designed to minimize contamination risks
- **Enhanced features:** These powder-free cleanroom gloves are manufactured with a beaded cuff, enhancing their strength and facilitating double-donning
- **Heightened tactility and comfort:** BioClean™ Ultimate BUPS polychloroprene gloves are also comfortable, with good tactility levels for engaging in precision work



### KEY FEATURES & BENEFITS

- **ASTM and EN standards:** Chemical and chemotherapy drug resistance
- **ESD properties:** Component integrity protection
- **Non-particulating EasyTear packaging:** Fewer contamination risks

### Industries

- Controlled and Critical Environments
- Production and Manufacturing
- Lab and R&D
- Sterile Processing Department

[WWW.ANSELL.COM](http://WWW.ANSELL.COM)





# Ultimate BUPS

Sterile Disposable Neoprene  
(Polychloroprene) Cleanroom Glove

## TECHNICAL DATA SHEET

Product Information	
Material	Neoprene (Polychloroprene)
Color	Natural
Shape	Hand Specific
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 14001, Manufacturing QMS Audit Standards ISO 9001, PPE Regulation 2016 425 Module D, NEBB Certified Cleanrooms
Regulatory/Standard Compliance	ASTM D6978, CE 0598, EN ISO 21420:2020, EN 421:2010, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 455 Part 2, Category III
Packaging	One pair per inner PE wallet; one wallet per sealed EasyTear PE pouch; 10 pouches per sealed outer PE bag; 20 outer bags per lined carton (200 pairs)
Storage	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country of Origin	Indonesia
Available sizes	6, 6.5, 7, 7.5, 8, 8.5, 9
Powder Content	Powder-Free
External Glove Surface	Textured Fingers and Palm
Internal Glove Surface	Polymer Coated
Sterilization Method	GAMMA irradiation (25 kGy)
Sterilization Minimum Dose	25kGy
Sterility Assurance Level	10 <sup>-6</sup>
Cleanroom Class	Class 10/ISO Class 4 & EU GMP Grade A/B and other sterile cleanrooms
Shelf Life	Five (5) years from date of manufacture.
Tested for use with Chemotherapy Drugs	Yes
Protein Level	N/A: contains no natural rubber latex
Anti-static	Yes



# Ultimate BUPS

Sterile Disposable Neoprene  
(Polychloroprene) Cleanroom Glove

Physical Properties								Testing Method
Sizes	6	6.5	7	7.5	8	8.5	9	
Typical Length (mm/in)	300 / 12							EN 420
Palm Width (mm/in)	77/3	83/3.2	89/3.5	95/3.7	102/4	108/4.2	115/4.5	
Freedom from Holes	0.65 AQL Performance Level 3							EN 374-2
Typical Particle Count ≥0.5µm (counts / cm <sup>2</sup> )	<1300							IEST-RP-CC005.4
Target Single Wall Palm Thickness (mm/mil)	0.11 / 4.33							EN 455-2
Target Single Wall Finger Thickness (mm/mil)	0.14 / 5.51							EN 455-2
Target Single Wall Cuff Thickness (mm/mil)	0.09 / 3.54							EN 455-2
Force at Break (N) Before Aging	≥ 9 N							EN 455-2
Force at Break (N) During Aging	≥ 9 N							EN 455-2

## IONIC CONTENT

Concentration in µg/cm <sup>2</sup>	Typical	Concentration in µg/cm <sup>2</sup>	Typical
Ammonium	0.005	Nitrate	0.795
Bromide	Not Detected	Nitrite	Not Detected
Calcium	0.582	Phosphate	Not Detected
Chloride	0.494	Potassium	0.378
Fluoride	Not Detected	Sodium	0.123
Lithium	Not Detected	Sulphate	0.041
Magnesium	0.025	Zinc	Not Tested

## ORDERING INFORMATION

SIZE	6	6.5	7	7.5	8	8.5	9
REORDER NO.	BUPS-60	BUPS-65	BUPS-70	BUPS-75	BUPS-80	BUPS-85	BUPS-90

## PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
T: +32 (0) 2 528 74 00  
F: +32 (0) 2 528 74 01

### Asia Pacific Region

Ansell Global Trading Center  
T: +603 8310 6688  
F: +603 8310 6699

### North America Region

Ansell Healthcare Products LLC  
US T: +1 800 800 0444  
US F: +1 800 800 0445  
CA T: +1-800-363-8340

### Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.  
T: +52 442 248 1544 / 248 3133

### Australia

Ansell Limited  
T: +61 1800 337 041  
F: +61 1800 803 578

### UK

Ansell Nitritex  
T: +44 1638 663338  
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:  
[www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

