

# GAMMEX<sup>®</sup> Non-Latex PI

S U R G I C A L   G L O V E S

## Soft, skin-friendly polyisoprene surgical gloves, for optimal comfort and allergy protection

- **Optimal allergy protection:** Being non-latex, these polyisoprene surgical gloves eliminate the risk of a latex allergy while minimizing sensitization to a Type IV chemical allergy as chemical accelerators known to cause skin allergies have been eliminated from its skin-friendly formulation
- **Superior comfort and dexterity:** Ansell's proprietary soft yet strong formulation delivers exceptional fit and feel
- **Excellent grip and sensitivity:** The smooth surface of these surgical gloves has a micro-textured finish to ensure confident handling and instrument grip during the most rigorous procedures
- **Effortless donning and doffing:** A unique internal coating aids dry and damp hand donning, either when first donned or during intra-operative changes
- **Chemo-tested:** Tested to ASTM D6978, these GAMMEX<sup>®</sup> sterile gloves are safe to use with a wide range of chemotherapy drugs\*\*



\*\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

### Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Extra comfort during long, sustained procedures
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies
- Having passed the Mouse Embryo Assay (MEA) test, may be used for In Vitro Fertilization (IVF) procedures



### FEATURED TECHNOLOGIES



## GAMMEX<sup>®</sup> Innovation that Protects

Driven by innovation, GAMMEX<sup>®</sup> surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Non-Latex PI

## SURGICAL GLOVES

### PRODUCT INFORMATION

|  |   |
|--|---|
| <b>Material</b>                                | PI (Polyisoprene)   |
| <b>Powder Content</b>                          | Powder-Free   |
| <b>Color</b>                                   | White   |
| <b>Thickness</b>                               | Standard  |
| <b>Cuff Style</b>                              | Beaded with SUREFIT™ Technology   |
| <b>External Glove Surface</b>                  | Smooth with Micro-Textured Finish   |
| <b>Internal Glove Surface</b>                  | Polymer Coated with DERMASHIELD™ Technology   |
| <b>Allergy Prevention</b>                      | Latex (Type I)  |
| <b>Grip Level</b>                              | Moderate  |
| <b>Double Gloving Recommendation</b>           | Outer glove or underglove   |
| <b>Tested For Use With Chemotherapy Drugs</b>  | Yes, in accordance with ASTM D6978 and US FDA cleared   |
| <b>Freedom from Holes (Inspection level I)</b> | 0.65 AQL  |
| <b>Sterile</b>                                 | Yes   |
| <b>Product Standards</b>                       | AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F739, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107 |
| <b>Quality / Environmental Standards</b>       | EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001  |
| <b>Packaging</b>                               | 50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case  |
| <b>Shelf Life</b>                              | 3 years   |
| <b>Storage Instructions</b>                    | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.  |
| <b>Vulcanization Chemical Accelerators</b>     | Xanthates, Zinc Dialkyldithiocarbamate  |

### ORDERING INFORMATION

| Size                | 5.5       | 6         | 6.5       | 7         | 7.5       | 8         | 8.5       | 9         |
|---------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| <b>Product Code</b> | 340056055 | 340056060 | 340056065 | 340056070 | 340056075 | 340056080 | 340056085 | 340056090 |

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: csjac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com