

Validated sterile nitrile RABS/isolator sleeves, ensuring low particle levels and excellent ESD properties

Port size: 152mm/6"

Thickness: 0.50mm/20 mil

- **Assured quality:** The BioClean™ Sterile Nitrile RABS/Isolator Sleeve GSL15NITPP26 is triple bagged, for assured integrity, and 100% visually inspected, ensuring high levels of quality
- **Enhanced versatility:** These nitrile protective sleeves are compatible with a wide range of cuff systems, extending protective coverage to the arm
- **Specialized protection:** In addition, they are autoclave compatible, offering wearers' arm protection in highly pressurized, high-temperature conditions
- **Added benefits:** These arm protection sleeves also feature extremely low levels of particles and offer excellent ESD properties



KEY FEATURES & BENEFITS

- **Triple bagged and 100% visually inspected:** Safety and quality assurance
- **Compatibility with varied cuff systems:** Truly versatile personal protection
- **Autoclave compatibility:** For pressurized, high-temperature conditions

Industries

- Pharmaceutical Manufacturing
- Laboratory Research and Development
- Raw Material and Preparation
- Biotechnology Manufacturing
- Sterile Processing Department
- Semiconductor manufacturing



GSL15NITPP26

Validated Sterile Nitrile RABS / Isolator Sleeve

TECHNICAL DATA SHEET

Product Reference	GSL15NITPP26
Material	Nitrile
Color	White
External Glove Surface	Smooth
Cuff Style	Beaded
Audit Standards	ISO 14001, ISO 9001:2008, NEBB Certified Cleanrooms
Packaging Overview	One piece per sealed inner PE bag; one inner bag per sealed second inner PE bag; one second inner bag per sealed outer PE bag; 20 outer bags per lined inner white Correx polyethylene box (20 pieces).
Storage Instructions	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.
Country Of Origin	Malaysia
Temperature Range	From -10°C to 120°C
Sterilization Method	GAMMA irradiation (25 kGy)
Cleanroom Class	Class 10/ISO 4 & EU GMP Grade A
Shelf Life	Three (3) years from date of manufacture.
Protein Level	N/A: contains no natural rubber latex
Antistatic	Yes
Standards	Category III, CE 0493, EC 1935/2004, ISO 11137-1:2006, EN 388:2016 +A1:2018, EN 420:2003 + A1:2009, EN 421:2010, EN 455 Part 2, EN ISO 374-1:2016, EN ISO 374-5:2016



GSL15NITPP26

Validated Sterile Nitrile RABS / Isolator Sleeve

PHYSICAL PROPERTIES

	Typical Values	Testing Method
Typical Length (mm/in)	660/26	EN ISO 21420
Port Size	6"	
Target Single Wall Cuff Thickness (mm/mil)	0.50/20	EN 455-2

ORDERING INFORMATION

	PORT SIZE	6"/152mm
	SIZE	6"/152mm
GSL15NITPP26	REORDER NO.	GSL15NITPP26

Technology



NO RUBBER LATEX

Performance Standards and Regulatory Compliance



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

