

## A revolution in skin moisturization with sensational grip

- **Skin moisturization:** MICROFLEX™ NEOSOFT™ are comfortable neoprene medical exam gloves designed with HYDRASOFT™, a unique technology that helps skin retain moisture to improve skin health. In fact, lab tests show wearing MICROFLEX™ NEOSOFT™ gloves help skin retain 2X more moisture\*
- **Elevated comfort:** The gloves' soft and stretchy formulation provides exceptional comfort and ease of use
- **Exceptional grip:** Advanced neoprene formulation and textured fingers deliver a secure, reliable wet grip
- **Tactile sensitivity:** An ultra-lightweight design allows for enhanced tactile sensitivity when handling delicate instruments, parts and pieces



\*based on internal testing (data on file)

### Industries

- Healthcare
- Dental

### Recommended For

- Enhanced grip in wet environments
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Extra comfort during long, sustained procedures
- Damp donning
- Standard patient examinations
- Protection from Type I latex allergy in HCW's or patients
- Helping protect skin from drying or cracking

# TECHNICAL DATA SHEET

Product Information	
<b>Material</b>	Neoprene
<b>External Glove Surface</b>	Textured Fingers
<b>Audit Standards</b>	ISO 13485, ISO 9001
<b>Product Certification</b>	ASTM D6977, EN ISO 21420:2020, FDA21 CFR 177-2600-US Food Contact Approved, US QSR/FDA 510(k) Medical Examination Grade, ISO 13485, ISO 9001
<b>Packaging Overview</b>	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
<b>Storage Instructions</b>	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
<b>Country Of Origin</b>	Indonesia
<b>Product Segmentation</b>	Standard Thin
<b>Antistatic</b>	No
<b>Silicone Free</b>	No

## PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	240 / 9.45		ASTM D3767,EN 420
Freedom from Holes (Inspection level I)	1.5 AQL		ASTM D5151,EN 374-2
Palm Thickness (mm/mils)	0.07 / 2.8		ASTM D3767,EN 420
Finger Thickness (mm/mils)	0.09 / 3.5		ASTM D3767,EN 420
	<b>BEFORE AGING</b>	<b>AFTER AGING</b>	
Elongation at Break (%)	≥ 500	≥ 400	ASTM D412

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	73837060	73837070	73837080	73837090	73837100

### Performance Standards and Regulatory Compliance



ISO 13485 ISO 9001

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

#### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV

T: +32 (0) 2 528 74 00

F: +32 (0) 2 528 74 01

#### Asia Pacific Region

MY: +603 8310 6688

apac.medical@ansell.com

CN: +86 21 38275000

infochina@ansell.com

JP: +813 5549 8151

info.medical.jp@ansell.com

#### North America Region

Ansell Healthcare Products LLC

T: +1 800 800 0444

F: +1 800 800 0445

#### Latin America & Caribbean Region

BR: +55 (11) 3356-3100

latam.medical@ansell.com

MX: +52 442 296 2050

latam.medical@ansell.com

CO: +57 1 288 3247

latam.medical@ansell.com

#### Australia

Ansell Limited

T: +61 1800 337 041

F: +61 1800 803 578

#### Russia

Ansell PYC

Ten. +7 495 258 13 16

Technology



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2025 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

