


Sterile Disposable Cleanroom Gloves, Extended Coverage

- **Heightened protection:** BioClean™ 73-425 sterile disposable cleanroom gloves feature a distinctive neoprene formulation that provides enhanced protection. They offer extended forearm coverage, and deliver high resistance to chemicals.
- **Enhanced Comfort:** These cleanroom gloves are designed with textured fingers and palms, made from an exceptionally soft material that ensures superior comfort, fit, feel, and flexibility.
- **Reduced Allergy risks:** These neoprene gloves contain no latex proteins and are chemical accelerator-free, preventing the risk of Type I and Type IV allergies and sensitivities.



KEY FEATURES & BENEFITS

- **Dermashield Technology:** Featuring proprietary inner coating for fast and easy damp or dry hand donning
- **Chemical accelerator free:** Type I and Type IV allergy risk eliminated
- **Longer length:** 400mm/16" length for extended protection of the forearm and elimination risk of skin exposure and gaps between glove and garment.

Industries

- Pharmaceutical Manufacturing
- Medical Device Manufacturing
- Laboratory Research and Development
- Biotechnology Manufacturing
- Sterile Processing Department

Recommended For

- Blending, compounding solids and liquids
- Packaging and labeling
- Transferring liquids and solids
- Biotechnologies Blending operations
- Preparing pharmaceutical products

TECHNICAL DATA SHEET

| Product Information | |
|--|--|
| Material | Neoprene (Polychloroprene) |
| Color | Brown |
| Shape | Hand Specific |
| Cuff | Beaded |
| Manufacturing/QMS Audit Standards | ISO 9001:2015, NEBB Certified Cleanrooms |
| Regulatory/Standard Compliance | CE 0493, EN 421:2010, EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016, Category III, UKCA |
| Packaging | <p>One pair glove per PE inner wallet, One wallet per sealed PE pouch, 25 pouches per PE inner bag, 2 inner bag per PE outer bag, 4 outer bag per carton with liner, 200 pairs per carton</p> <p>More sustainable packaging: These sterile gloves are packed in recyclable plastic packaging* and delivered in recycled cardboard shipper cases.</p> <p>*Inner wrap, pouch, bag and liner are made from polyethylene (PE) based film. Always check your local recyclable status as these materials may not be considered suitable for recycling in your location.</p> |
| Storage | Keep away from direct sunlight; store in a dry place and keep in the original packaging. Keep away from ozone sources. If products are properly stored, as indicated, they won't lose their performances or change characteristics significantly. If products could be affected by ageing or storage, the expiry date is mentioned on the packaging materials. |
| Country of Origin | Malaysia |
| Available sizes | 6, 6.5, 7, 7.5, 8, 8.5, 9, 10 |
| Powder Content | Powder-Free |
| External Glove Surface | Textured Fingers and Palm |
| Internal Glove Surface | Polymer Coated with DERMASHIELD™ Technology |
| Sterilization Method | GAMMA irradiation (25 kGy) |
| Sterilization Minimum Dose | 25kGy |
| Sterility Assurance Level | 10 ⁻⁶ |
| Cleanroom Class | Class 100/ISO Class 5 & EU GMP Grade A/B and other sterile cleanrooms |
| Shelf Life | Three (3) years from date of manufacture. |
| Tested for use with Chemotherapy Drugs | No |
| Protein Level | N/A: contains no natural rubber latex |
| Anti-static | No |



73-425

Sterile Disposable Neoprene (Polychloroprene) Cleanroom Glove

| Physical Properties | Sizes | | | | | | | | Testing Method |
|--|-----------|--------|--------|--------|---------|---------|---------|---------|-----------------|
| | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 | 10 | |
| Typical Length (mm/in) | 400 /16 | | | | | | | | EN ISO 21420 |
| Palm Width (mm/in) | 81/3.2 | 88/3.5 | 95/3.7 | 99/3.9 | 106/4.2 | 111/4.4 | 117/4.6 | 123/4.8 | |
| Freedom from Holes | 1.5 AQL | | | | | | | | EN 374-2 |
| Typical Particle Count $\geq 0.5\mu\text{m}$ (counts / cm ²) | <3500 | | | | | | | | TEST-RP-CC005.4 |
| Target Single Wall Palm Thickness (mm/mil) | 0.20/7.87 | | | | | | | | EN 455-2 |
| Target Single Wall Finger Thickness (mm/mil) | 0.22/8.66 | | | | | | | | EN 455-2 |
| Target Single Wall Cuff Thickness (mm/mil) | 0.17/6.69 | | | | | | | | EN 455-2 |
| Ultimate tensile strength (MPa) During Aging | ≥ 14 | | | | | | | | ASTM D412-06a |
| Force at Break (N) During Aging | ≥ 9 | | | | | | | | EN 455-2 |

ORDERING INFORMATION

| SIZE | 6 | 6.5 | 7 | 7.5 | 8 | 8.5 | 9 | 10 |
|-------------|----------|----------|----------|----------|----------|----------|----------|----------|
| REORDER NO. | 73425060 | 73425065 | 73425070 | 73425075 | 73425080 | 73425085 | 73425090 | 73425100 |

PERFORMANCE STANDARDS AND REGULATORY COMPLIANCE



For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
US T: +1 800 800 0444
US F: +1 800 800 0445
CA T: +1-800-363-8340

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

UK

Ansell Nitritex
T: +44 1638 663338
F: +44 1638 668890

Technology



Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

